

## ADDING THE SUPER TO HEROES

Ever wondered what makes a superhero? What makes you act by heart and without fear? Or react on instinct without thought of consequence?

We believe it's certainty. Knowing that when you are fully geared-up, you become 100% of what you can be. Knowing that whatever you put your heart into, you will do it with maximized power. Knowing that the result you get today will be the best – the very best. At least for today.

Rehband's Compression products are said to be the best compression gear, outclassing the competition in all areas. Using it won't give you the ability to fly or grow the strength to lift a train. But with Rehband, you can be certain you'll be all that you can be.



"On my mountain bike, I need to be sure that the equipment I use will aid my training, racing, recovery and trailblazing activities and commitments. It also needs to last as long as I do on my rides in some of the toughest environments known to cycling.

This is why I use Rehband Compression Garments. They keep my body feeling fresh, smelling fresh and cycling for long periods without aches or pains. Even lactate accumulates much later.

Without Rehband, I seriously doubt if I could have completed my 59hr and 36minute non-stop mountain bike ride in which I covered 721km and climbed 22300m through Wales' toughest countryside.

– David Buchanan  
Guinness World Record™ holder: "The greatest distance achieved on a mountain bike in 48hrs"



## WHAT REALLY MATTERS TO US IS WHAT YOU THINK.

All Rehband products are tested by sportsmen for many different types of sports. We want to know that our products can be used for a marathon, for a footballer, athletics, under special conditions for martial arts, under extreme conditions as David Buchanan's World Record attempt as well as being used for indoor sports as basket ball and handball, and many more sports. It's confirmed from all these tests that we have a superior product compared to our

competitors. The response that we got was telling us that we were beating all competition in all areas and especially the most important ones as compression, moisture control and shock absorption. The combination of all long term testing with the fact that it's tested on so many different sports and also the wide range of athletes, from elite to amateurs, gives you products that is tested as The no 1 Compression Gear.




"I have been working with professional athletes and national teams for years and realize that the recuperation between matches or physical endeavours, especially the first 24 hours, makes the difference in performance and injury risk. I have tried many methods, however, compression has always been an integral part and the difficulty is getting good compression that shapes with the body without being too hot or tight. Compression socks have been widely used but the rest of the body has been difficult. We have had to rely on different submersion baths which are not effective in the long run and different types of circulation training. When the Rehband Compression products were introduced, the players felt that they were actively getting ready for the next endeavor immediately after a match."

– Dale Reese  
Physiotherapist

## ✓ FIRST-CLASS COMPRESSION

The results from Rehband's tests show higher compression and better sustainability thanks to more technologically advanced materials



## ✓ OUTSTANDING ANATOMICAL FIT

Our garments offer tight fitting, high performance fabric, uniquely composed to allow the user to be comfortable and maintain the agility and mobility needed during tenuous exercise




## ✓ UNBEATEN MOISTURE CONTROL

Staying dry whilst exercising or competing is essential for maximum performance. Our compression garments' unique composition channel the moisture away from the body, keeping you comfortable.



## ✓ UNBEATABLE SHOCK ABSORPTION

Rehband's Protection Gear with the unique, larger total pelot area differs from other protection on the market. The Rehband pelot spreads the total impact force evenly into many small impact forces, instead of unevenly spreading the impact via seperated pelot elements, as is often the case with the commonly used individual foam pad pelot.



REHband